How to use bimaps from the ".db" annotation packages

Marc Carlson, Herve Pages, Seth Falcon, Nianhua Li

October 1, 2014

1 Introduction

1.0.1 Purpose

AnnotationDbi is used primarily to create mapping objects that allow easy access from R to underlying annotation databases. As such, it acts as the R interface for all the standard annotation packages. Underlying each AnnotationDbi supported annotation package is at least one (and often two) annotation databases. AnnotationDbi also provides schemas for theses databases. For each supported model organism, a standard gene centric database is maintained from public sources and is packaged up as an appropriate organism or "org" package.

1.0.2 Database Schemas

For developers, a lot of the benefits of having the information loaded into a real database will require some knowledge about the database schema. For this reason the schemas that were used in the creation of each database type are included in AnnotationDbi. The currently supported schemas are listed in the DBschemas directory of AnnotationDbi. But it is also possible to simply print out the schema that a package is currently using by using its "_dbschema" method.

There is one schema/database in each kind of package. These schemas specify which tables and indices will be present for each package of that type. The schema that a particular package is using is also listed when you type the name of the package as a function to obtain quality control information.

The code to make most kinds of the new database packages is also included in AnnotationDbi. Please see the vignette on SQLForge for more details on how to make additional database packages.

1.0.3 Internal schema Design of org packages

The current design of the organism packages is deliberately simple and gene centric. Each table in the database contains a unique kind of information and also an internal identifier called _id. The internal _id has no meaning outside of the context of a single database. But _id does connect all the data within a single database.

As an example if we wanted to connect the values in the genes table with the values in the kegg table, we could simply join the two tables using the internal _id column. It is very important to note however that _id does not have any absolute significance. That is, it has no meaning outside of the context of the database where it is used. It is tempting to think that an _id could have such significance because within a single database,

it looks and behaves similarly to an entrez gene ID. But _id is definitely NOT an entrez gene ID. The entrez gene IDs are in another table entirely, and can be connected to using the internal _id just like all the other meaningful information inside these databases. Each organism package is centered around one type of gene identifier. This identifier is found as the gene_id field in the genes table and is both the central ID for the database as well as the foreign key that chip packages should join to.

The chip packages are 'lightweight', and only contain information about the basic probe to gene mapping. You might wonder how such packages can provide access to all the other information that they do. This is possible because all the other data provided by chip packages comes from joins that are performed by AnnotationDbi behind the scenes at run time. All chip packages have a dependency on at least one organism package. The name of the organism package being depended on can be found by looking at its "ORGPKG" value. To learn about the schema from the appropriate organism package, you will need to look at the "_dbschema" method for that package. In the case of the chip packages, the gene_id that in these packages is mapped to the probe_ids, is used as a foreign key to the appropriate organism package.

Specialized packages like the packages for GO and KEGG, will have their own schemas but will also adhere to the use of an internal _id for joins between their tables. As with the organism packages, this _id is not suitable for use as a foreign key.

For a complete listing of the different schemas used by various packages, users can use the available.dbschemas function. This list will also tell you which model organisms are supported.

```
require(org.Hs.eg.db)
## Loading required package: org.Hs.eg.db
## Loading required package:
                              {\it Annotation Dbi}
## Loading required package: BiocGenerics
## Loading required package:
                              parallel
##
## Attaching package: 'BiocGenerics'
##
## The following objects are masked from 'package:parallel':
##
##
      clusterApply, clusterApplyLB, clusterCall, clusterEvalQ, clusterExport,
      clusterMap, parApply, parCapply, parLapply, parLapplyLB, parRapply,
##
##
      parSapply, parSapplyLB
##
## The following object is masked from 'package:stats':
##
##
      xtabs
##
## The following objects are masked from 'package:base':
##
##
      anyDuplicated, append, as.data.frame, as.vector, cbind, colnames, do.call,
      duplicated, eval, evalq, Filter, Find, get, intersect, is.unsorted, lapply,
##
##
      Map, mapply, match, mget, order, paste, pmax, pmax.int, pmin, pmin.int,
      Position, rank, rbind, Reduce, rep. int, rownames, sapply, setdiff, sort,
##
##
      table, tapply, union, unique, unlist
##
## Loading required package: Biobase
```

```
## Welcome to Bioconductor
##
##
      Vignettes contain introductory material; view with 'browseVignettes()'.
      cite Bioconductor, see 'citation("Biobase")', and for packages
##
      'citation("pkqname")'.
##
##
## Loading required package:
                              GenomeInfoDb
## Loading required package:
                              DBI
require(AnnotationForge)
## Loading required package: AnnotationForge
available.dbschemas()
```

2 Examples

2.0.4 Basic information

The AnnotationDbi package provides an interface to SQLite-based annotation packages. Each SQLite-based annotation package (identified by a ".db" suffix in the package name) contains a number of AnnDbBimap objects in place of the environment objects found in the old-style environment-based annotation packages. The API provided by AnnotationDbi allows you to treat the AnnDbBimap objects like environment instances. For example, the functions [[, get, mget, and ls all behave the same as they did with the older environment based annotation packages. In addition, new methods like [, toTable, subset and others provide some additional flexibility in accessing the annotation data.

```
library("hgu95av2.db")
##
```

The same basic set of objects is provided with the db packages:

```
ls("package:hgu95av2.db")
    [1] "hgu95av2"
                                 "hgu95av2_dbconn"
##
                                 "hgu95av2_dbInfo"
##
    [3] "hgu95av2_dbfile"
    [5] "hgu95av2_dbschema"
                                 "hgu95av2.db"
##
    [7] "hgu95av2ACCNUM"
                                 "hgu95av2ALIAS2PROBE"
##
    [9] "hgu95av2CHR"
                                 "hgu95av2CHRLENGTHS"
##
## [11] "hgu95av2CHRLOC"
                                 "hgu95av2CHRLOCEND"
## [13] "hgu95av2ENSEMBL"
                                 "hgu95av2ENSEMBL2PR0BE"
## [15] "hgu95av2ENTREZID"
                                 "hgu95av2ENZYME"
## [17] "hgu95av2ENZYME2PR0BE"
                                 "hgu95av2GENENAME"
## [19] "hgu95av2GO"
                                 "hgu95av2G02ALLPR0BES"
## [21] "hgu95av2G02PR0BE"
                                 "hgu95av2MAP"
## [23] "hgu95av2MAPCOUNTS"
                                 "hgu95av20MIM"
## [25] "hgu95av20RGANISM"
                                 "hgu95av20RGPKG"
## [27] "hgu95av2PATH"
                                 "hgu95av2PATH2PR0BE"
```

```
## [29] "hgu95av2PFAM" "hgu95av2PMID"

## [31] "hgu95av2PMID2PROBE" "hgu95av2PROSITE"

## [33] "hgu95av2REFSEQ" "hgu95av2SYMBOL"

## [35] "hgu95av2UNIGENE" "hgu95av2UNIPROT"
```

Exercise 1

Start an R session and use the library function to load the hgu95av2.db software package. Use search() to see that an organism package was also loaded and then use the approriate "_dbschema" methods to the schema for the hgu95av2.db and org.Hs.eg.db packages.

It is possible to call the package name as a function to get some QC information about it.

```
qcdata = capture.output(hgu95av2())
## hqu95av2PFAM is deprecated because up to date IPI IDs are no longer available.
## Please use select() if you need access to PFAM or PROSITE accessions.
##
## hqu95av2PROSITE is deprecated because up to date IPI IDs are no longer available.
## Please use select() if you need access to PFAM or PROSITE accessions.
head(qcdata, 20)
    [1] "Quality control information for hgu95av2:"
    [2] ""
##
  [3] ""
##
## [4] "This package has the following mappings:"
## [5] ""
    [6] "hgu95av2ACCNUM has 12625 mapped keys (of 12625 keys)"
    [7] "hgu95av2ALIAS2PROBE has 33560 mapped keys (of 103510 keys)"
## [8] "hgu95av2CHR has 11531 mapped keys (of 12625 keys)"
## [9] "hgu95av2CHRLENGTHS has 93 mapped keys (of 93 keys)"
## [10] "hgu95av2CHRLOC has 11466 mapped keys (of 12625 keys)"
## [11] "hgu95av2CHRLOCEND has 11466 mapped keys (of 12625 keys)"
## [12] "hgu95av2ENSEMBL has 11436 mapped keys (of 12625 keys)"
## [13] "hgu95av2ENSEMBL2PROBE has 9756 mapped keys (of 28046 keys)"
## [14] "hgu95av2ENTREZID has 11533 mapped keys (of 12625 keys)"
## [15] "hgu95av2ENZYME has 2125 mapped keys (of 12625 keys)"
## [16] "hgu95av2ENZYME2PROBE has 786 mapped keys (of 975 keys)"
## [17] "hgu95av2GENENAME has 11533 mapped keys (of 12625 keys)"
## [18] "hgu95av2GO has 11228 mapped keys (of 12625 keys)"
## [19] "hgu95av2G02ALLPROBES has 16526 mapped keys (of 18078 keys)"
## [20] "hgu95av2G02PR0BE has 12422 mapped keys (of 14134 keys)"
```

Alternatively, you can get similar information on how many items are in each of the provided maps by looking at the MAPCOUNTs:

```
hgu95av2MAPCOUNTS
```

To demonstrate the environment API, we'll start with a random sample of probe set IDs.

```
all_probes <- ls(hgu95av2ENTREZID)
length(all_probes)
## [1] 12625
set.seed(0xa1beef)
probes <- sample(all_probes, 5)
probes
## [1] "31882_at" "38780_at" "37033_s_at" "1702_at" "31610_at"</pre>
```

The usual ways of accessing annotation data are also available.

```
hgu95av2ENTREZID[[probes[1]]]

## [1] "9136"

hgu95av2ENTREZID$"31882_at"

## [1] "9136"

syms <- unlist(mget(probes, hgu95av2SYMBOL)))

syms

## 31882_at 38780_at 37033_s_at 1702_at 31610_at

## "RRP9" "AKR1A1" "GPX1" "IL2RA" "PDZK1IP1"
```

The annotation packages provide a huge variety of information in each package. Some common types of information include gene symbols (SYMBOL), GO terms (GO), KEGG pathway IDs (KEGG), ENSEMBL IDs (ENSEMBL) and chromosome start and stop locations (CHRLOC and CHRLOCEND). Each mapping will have a manual page that you can read to describe the data in the mapping and where it came from.

?hgu95av2CHRLOC

Exercise 2

For the probes in 'probes' above, use the annotation mappings to find the chromosome start locations.

2.0.5 Manipulating Bimap Objects

Many filtering operations on the annotation *Bimap* objects require conversion of the *AnnDbBimap* into a *list*. In general, converting to lists will not be the most efficient way to filter the annotation data when using a SQLite-based package. Compare the following two examples for how you could get the 1st ten elements of the hgu95av2SYMBOL mapping. In the 1st case we have to get the entire mapping into list form, but in the second case we first subset the mapping object itself and this allows us to only convert the ten elements that we care about.

```
system.time(as.list(hgu95av2SYMBOL)[1:10])
## vs:
system.time(as.list(hgu95av2SYMBOL[1:10]))
```

There are many different kinds of *Bimap* objects in AnnotationDbi, but most of them are of class *AnnDbBimap*. All /RclassBimap objects represent data as a set of left and right keys. The typical usage of these mappings

is to search for right keys that match a set of left keys that have been supplied by the user. But sometimes it is also convenient to go in the opposite direction.

The annotation packages provide many reverse maps as objects in the package name space for backwards compatibility, but the reverse mappings of almost any map is also available using revmap. Since the data are stored as tables, no extra disk space is needed to provide reverse mappings.

So now that you know about the revmap function you might try something like this:

```
as.list(revmap(hgu95av2PATH)["00300"])
## $`00300`
## [1] "36132_at" "35870_at"
```

Note that in the case of the PATH map, we don't need to use revmap(x) because hgu95av2.db already provides the PATH2PROBE map:

```
x <- hgu95av2PATH
## except for the name, this is exactly revmap(x)
revx <- hgu95av2PATH2PROBE
revx2 <- revmap(x, objName="PATH2PROBE")
revx2
## PATH2PROBE map for chip hgu95av2 (object of class "ProbeAnnDbBimap")
identical(revx, revx2)
## [1] TRUE
as.list(revx["00300"])
## $`00300`
## [1] "36132_at" "35870_at"</pre>
```

Note that most maps are reversible with revmap, but some (such as the more complex GO mappings), are not. Why is this? Because to reverse a mapping means that there has to be a "value" that will always become the "key" on the newly reversed map. And GO mappings have several distinct possibilities to choose from (GO ID, Evidence code or Ontology). In non-reversible cases like this, AnnotationDbi will usually provide a pre-defined reverse map. That way, you will always know what you are getting when you call revmap

While we are on the subject of GO and GO mappings, there are a series of special methods for GO mappings that can be called to find out details about these IDs. Term,GOID, Ontology, Definition,Synonym, and Secondary are all useful ways of getting additional information about a particular GO ID. For example:

```
Term("GD:0000018")
## Loading required package: GD.db
## GD:0000018
## "regulation of DNA recombination"
Definition("GD:0000018")
```

```
##
## "Any process that modulates the frequency, rate or extent of DNA recombination, a DNA metaboli
```

Exercise 3

Given the following set of RefSeq IDs: c("NG_005114","NG_007432","NG_008063"), Find the Entrez Gene IDs that would correspond to those. Then find the GO terms that are associated with those entrez gene IDs. org.Hs.eg.db packages.

2.0.6 The Contents and Structure of Bimap Objects

Sometimes you may want to display or subset elements from an individual map. A *Bimap* interface is available to access the data in table (*data.frame*) format using [and toTable.

```
head(toTable(hgu95av2G0[probes]))
##
     probe_id
                   go_id Evidence Ontology
## 1 1702_at GD:0006915
                               TAS
                                         BP
## 2 1702_at GO:0006955
                               TAS
                                         BP
## 3 1702_at GO:0007166
                               TAS
                                         BP
## 4 1702_at GO:0008283
                               TAS
                                         BP
## 5 1702_at GO:0002437
                                         BP
                               IEA
## 6 1702_at GO:0006924
                               IEA
                                         BP
```

The toTable function will display all of the information in a *Bimap*. This includes both the left and right values along with any other attributes that might be attached to those values. The left and right keys of the *Bimap* can be extracted using Lkeys and Rkeys. If is is necessary to only display information that is directly associated with the left to right links in a *Bimap*, then the links function can be used. The links returns a data frame with one row for each link in the bimap that it is applied to. It only reports the left and right keys along with any attributes that are attached to the edge between these two values.

Note that the order of the cols returned by toTable does not depend on the direction of the map. We refer to it as an 'undirected method':

```
toTable(x)[1:6, ]
     probe_id path_id
## 1 38187_at
                00232
## 2 38187_at
                00983
## 3 38187_at
                01100
## 4 38912_at
                00232
## 5 38912_at
                00983
## 6 38912_at
                01100
toTable(revx)[1:6, ]
     probe_id path_id
## 1 38187_at
                00232
## 2 38187_at
                00983
## 3 38187_at
                01100
                00232
## 4 38912_at
```

```
## 5 38912_at 00983
## 6 38912_at 01100
```

Notice however that the Lkeys are always on the left (1st col), the Rkeys always in the 2nd col

There can be more than 2 columns in the returned data frame:

3 cols:

```
toTable(hgu95av2PFAM)[1:6, ] # the right values are tagged

## probe_id ipi_id PfamId

## 1 1000_at IPI00018195 PF00069

## 2 1000_at IPI00984821 PF00069

## 3 1001_at IPI00019530 PF00041

## 4 1001_at IPI00019530 PF07714

## 5 1001_at IPI00019530 PF12661

## 6 1001_at IPI01009976 PF00041

as.list(hgu95av2PFAM["1000_at"])

## $`1000_at`

## IPI00018195 IPI00984821

## "PF00069" "PF00069"
```

But the Rkeys are ALWAYS in the 2nd col.

For length() and keys(), the result does depend on the direction, hence we refer to these as 'directed methods':

```
length(x)
## [1] 12625
length(revx)
## [1] 229
allProbeSetIds <- keys(x)
allKEGGIds <- keys(revx)</pre>
```

There are more 'undirected' methods listed below:

Notice how they give the same result for x and revmap(x)

You might be tempted to think that Lkeys and Llength will tell you all that you want to know about the

left keys. But things are more complex than this, because not all keys are mapped. Often, you will only want to know about the keys that are mapped (ie. the ones that have a corresponding Rkey). To learn this you want to use the mappedkeys or the undirected variants mappedkeys and mappedkeys. Similarly, the count.mappedkeys, count.mappedkeys and count.mappedkeys methods are very fast ways to determine how many keys are mapped. Accessing keys like this is usually very fast and so it can be a decent strategy to subset the mapping by 1st using the mapped keys that you want to find.

```
x = hgu95av2ENTREZID[1:10]
## Directed methods
mappedkeys(x)
                       # mapped keys
## [1] "1000_at"
                  "1001_at" "1002_f_at" "1003_s_at" "1004_at"
## [6] "1005_at"
                  "1006_at" "1008_f_at" "1009_at"
count.mappedkeys(x)
                    # nb of mapped keys
## [1] 9
## Undirected methods
mappedLkeys(x)
                       # mapped left keys
## [1] "1000_at"
                              "1002_f_at" "1003_s_at" "1004_at"
                  "1001_at"
## [6] "1005_at"
                  "1006_at"
                             "1008_f_at" "1009_at"
count.mappedLkeys(x) # nb of mapped Lkeys
## [1] 9
```

If you want to find keys that are not mapped to anything, you might want to use isNA.

```
y = hgu95av2ENTREZID[isNA(hgu95av2ENTREZID)]  # usage like is.na()
Lkeys(y)[1:4]
## [1] "1007_s_at" "1047_s_at" "1089_i_at" "108_g_at"
```

Exercise 4

How many probesets do not have a GO mapping for the hgu95av2.db package? How many have no mapping? Find a probeset that has a GO mapping. Now look at the GO mappings for this probeset in table form.

2.0.7 Some specific examples

Lets use what we have learned to get information about the probes that are are not assigned to a chromosome:

```
x <- hgu95av2CHR
Rkeys(x)
## [1] "19" "12" "8" "14" "3" "2" "17" "16" "9" "X" "6" "1" "7"
## [14] "10" "11" "22" "5" "18" "15" "Y" "20" "21" "4" "13" "MT" "Un"
chroms <- Rkeys(x)[23:24]
chroms
## [1] "4" "13"</pre>
```

```
Rkeys(x) <- chroms</pre>
toTable(x)
##
         probe_id chromosome
## 1
        1029_s_at
## 2
         1036_at
                           4
## 3
          1058_at
                           13
## 4
          1065_at
                          13
## 5
          1115_at
                           4
## 6
          1189_at
                          13
## 7
          1198_at
                          13
          1219_at
                           4
## 8
## 9
        1220_g_at
                           4
          1249_at
                            4
## 10
          1285_at
## 11
## 12
          1303_at
                           4
## 13
          1325_at
                           4
## 14
        1348_s_at
                           13
## 15
        1369_s_at
                           4
                           4
## 16
        1377_at
                           4
## 17
        1378_g_at
## 18
        1451_s_at
                          13
## 19
        1503_at
                           13
## 20
        1507_s_at
                           4
## 21
        1527_s_at
                          13
## 22
         1528_at
                          13
## 23
         1529_at
                          13
        1530_g_at
## 24
                           13
## 25
        1531_at
                           13
## 26
        1532_g_at
                           13
## 27
        1538_s_at
                           4
## 28
         1542_at
                           4
## 29
        1545_g_at
                           13
## 30
        1567_at
                           13
## 31
        1570_f_at
                           13
## 32
        1571_f_at
                          13
## 33
         1593_at
                           4
## 34
          1597_at
                          13
## 35
       1598_g_at
                           13
## 36
          159_at
                           4
## 37
          1600_at
                           4
## 38
         1604_at
                           4
## 39
                           4
        1605_g_at
## 40
          1616_at
                           13
## 41
                           4
          1624_at
## 42
        1629_s_at
                           4
## 43
        1670_at
                           13
        1672_f_at
## 44
                           13
```

##		1679_at	4
##	46	1708_at	4
##	47	1709_g_at	4
##	48	170_at	13
##	49	1720_at	4
##	50	1721_g_at	4
##	51	1731_at	4
##		1732_at	4
##		1819_at	13
##		1828_s_at	4
##		1836_at	4
##		1883_s_at	4
##		1888_s_at	4
##		1900_at	13
##		1905_s_at	13
##	60	1913_at	4
##	61	1914_at	13
##	62	1931_at	13
##	63	1934_s_at	4
##		1943_at	4
##		1954_at	4
##		1963_at	13
##		1964_g_at	13
##		1987_at	4
##		1988_at	4
##		1989_at	13
##		1990_g_at	13
##		2044_s_at	13
##		2062_at	4
##	74	2092_s_at	4
##	75	214_at	4
##	76	215_g_at	4
##		252_at	13
##		253_g_at	13
##		260_at	4
##		281_s_at	4
##		31314_at	4
##		31320_at	13
##		31333_at	4
##		31345_at	4
##		31349_at	4
##		31356_at	4
##	87	31382_f_at	4
##	88	31404_at	13
##	89	31408_at	4
##	90	31464_at	13
	91	31465_g_at	13
		Ü	

```
## 92 31516_f_at
                         13
      31543_at
                         4
## 93
## 94
      31562_at
                         13
## 95
      31584_at
                         13
## 96
      31628_at
                         13
## 97 31631_f_at
                         4
## 98 31639_f_at
                         13
## 99 31640_r_at
                         13
## 100 31670_s_at
                          4
                          4
## 101
        31684_at
## 102 31686_at
                          4
## 103
                          4
      31706_at
                          4
## 104 31744_at
## 105
      31753_at
                         13
## 106 31790_at
                         13
## 107
      31792_at
                          4
## 108 31805_at
                         4
## 109 31811_r_at
                         4
## 110
       31847_at
                         13
## 111
        31849_at
                         13
## 112
       31851_at
                         13
## 113 31876_r_at
                         4
                          4
## 114
        31894_at
## 115 31969_i_at
                          4
## 116 31970_r_at
                          4
## 117 32006_r_at
                          4
## 118 32026_s_at
                          4
## 119
        32080_at
                          4
## 120
      32102_at
                         13
## 121
        32145_at
                         4
## 122 32146_s_at
                         4
## 123
        32147_at
                         13
## 124
        32148_at
                         13
## 125 32163_f_at
                         4
                         4
## 126 32180_s_at
## 127
        32220_at
                         13
## 128
                          4
        32299_at
## 129 32349_at
## 130
      32353_at
                          4
                          4
## 131
      32357_at
## 132
        32368_at
                         13
## 133 32393_s_at
                         4
## 134 32439_at
                         13
        32446_at
                          4
## 135
## 136 32449_at
                          4
## 137
        32465_at
                         4
## 138
        32482_at
                         13
```

```
## 139
        32506_at
## 140
        32507_at
                           4
## 141
        32570_at
                           4
## 142
        32580_at
                           4
## 143
        32595_at
                           4
## 144
        32602_at
                           4
## 145
         32641_at
                          13
## 146
        32675_at
                           4
## 147
        32703_at
                           4
## 148
        32768_at
                          13
## 149
        32769_at
                           4
                           4
## 150
        32770_at
## 151
                           4
        32771_at
## 152
        32812_at
                           4
## 153
       32822_at
                           4
## 154
        32832_at
                           4
## 155
       32862_at
                          13
## 156
        32906_at
                          13
## 157
        32979_at
                           4
## 158 32986_s_at
                          13
## 159
        32998_at
                           4
## 160
        33013_at
                           4
                           4
## 161
         33050_at
## 162 33068_f_at
                           4
## 163 33069_f_at
                           4
## 164
         33100_at
                           4
## 165
                           4
         33150_at
## 166 33151_s_at
                           4
## 167
        33155_at
                           4
## 168
        33156_at
                           4
## 169
        33168_at
                          13
## 170 33171_s_at
                           4
## 171
         33172_at
                           4
## 172 33173_g_at
                           4
## 173
        33199_at
                          13
## 174
       33208_at
                          13
## 175
                           4
        33241_at
## 176
        33249_at
                           4
## 177
        33267_at
                           4
## 178
        33276_at
                          13
                           4
## 179
        33299_at
## 180
        33318_at
                          13
                           4
## 181
        33356_at
## 182
        33359_at
                           4
## 183
         33369_at
                           4
## 184 33370_r_at
                           4
## 185
                           4
         33382_at
```

```
## 186
         33483_at
## 187
                            4
         33488_at
## 188
         33490_at
                            4
## 189
         33494_at
                            4
## 190
         33519_at
                            4
## 191
         33520_at
                           13
## 192
         33525_at
                            4
## 193
         33526_at
                            4
## 194
         33529_at
                            4
## 195
         33536_at
                            4
## 196
         33544_at
                            4
## 197
                            4
         33564_at
## 198
         33576_at
                           13
## 199
         33584_at
                            4
## 200
         33596_at
                            4
## 201
         33657_at
                            4
## 202 33672_f_at
                            4
## 203 33673_r_at
                           4
## 204
         33687_at
                           13
## 205
         33700_at
                           13
## 206
         33733_at
                           4
## 207
        33791_at
                           13
## 208
        33823_at
                           4
## 209
        33827_at
                           13
## 210
                           4
        33837_at
## 211
        33859_at
                           13
## 212
         33975_at
                            4
## 213
         33990_at
                            4
## 214 33991_g_at
                            4
## 215
         33992_at
                            4
## 216
         33997_at
                            4
## 217
         34021_at
                            4
## 218
        34022_at
                            4
## 219
        34026_at
                           13
## 220
        34029_at
                           4
## 221
                           4
         34048_at
## 222
                           13
         34051_at
## 223
        34058_at
                            4
## 224
         34075_at
                            4
## 225
                            4
        34122_at
                            4
## 226
        34131_at
## 227
        34144_at
                            4
## 228
                            4
         34145_at
## 229
         34149_at
                            4
## 230 34170_s_at
                            4
## 231
         34181_at
                            4
## 232
                            4
         34198_at
```

```
## 233
        34211_at
                          13
## 234
                          13
         34239_at
## 235 34240_s_at
                          13
## 236
         34247_at
                           4
## 237
         34248_at
                           4
## 238 34275_s_at
                           4
## 239
        34284_at
                          13
## 240
       34307_at
                          13
## 241
       34319_at
                          4
## 242
       34324_at
                          13
## 243 34334_at
                          13
## 244
                          13
       34335_at
## 245
                           4
       34341_at
## 246 34342_s_at
                           4
## 247
        34353_at
                          4
## 248
        34398_at
                          13
## 249 34411_at
                          4
## 250
        34423_at
                          4
## 251
        34459_at
                          13
                           4
## 252 34476_r_at
                           4
## 253
        34482_at
## 254
        34512_at
                           4
                           4
## 255
       34551_at
## 256
       34564_at
                           4
## 257
       34565_at
                           4
## 258
       34578_at
                          13
## 259
        34583_at
                          13
## 260
                          4
        34596_at
## 261 34637_f_at
                          4
## 262 34638_r_at
                          4
## 263
        34657_at
                          13
## 264
        34672_at
                          13
## 265
       34745_at
                          4
## 266
         34803_at
                          13
## 267 34953_i_at
                          4
## 268 34954_r_at
                          4
## 269
                          13
        34955_at
## 270
        34973_at
                           4
## 271
        34984_at
                           4
## 272
                           4
       34988_at
                           4
## 273
        35020_at
## 274
        35021_at
                           4
## 275
                           4
        35025_at
## 276
        35028_at
                           4
## 277
        35039_at
                           4
## 278
        35053_at
                           4
## 279
                           4
        35061_at
```

```
## 280
         35063_at
                           4
## 281
         35081_at
                           13
## 282
         35105_at
                           13
## 283
         35107_at
                           13
## 284
         35110_at
                           13
## 285
         35131_at
                            4
## 286
         35134_at
                            4
## 287
         35140_at
                           13
## 288
         35147_at
                           13
## 289
                            4
         35164_at
## 290
         35181_at
                            4
## 291 35182_f_at
                            4
## 292
         35193_at
                           13
## 293
         35213_at
                           13
## 294
         35214_at
                            4
## 295
         35215_at
                            4
## 296
         35220_at
                            4
## 297
         35285_at
                            4
## 298
        35306_at
                            4
## 299
         35344_at
                           13
         35356_at
## 300
                            4
## 301
         35357_at
                            4
                            4
## 302
         35371_at
## 303 35372_r_at
                            4
## 304
         35400_at
                           13
## 305
         35410_at
                            4
## 306 35435_s_at
                            4
## 307
         35437_at
                            4
## 308
         35469_at
                           13
## 309
         35470_at
                           13
## 310 35471_g_at
                           13
## 311
         35481_at
                           13
                           4
## 312
         35507_at
## 313
         35523_at
                            4
## 314 35554_f_at
                           13
## 315 35555_r_at
                           13
## 316
         35591_at
                            4
## 317
         35656_at
                           13
## 318
        35662_at
                            4
## 319
                            4
        35664_at
                            4
## 320
        35678_at
## 321
        35689_at
                            4
## 322
                            4
        35698_at
## 323
         35725_at
                           13
## 324
         35730_at
                            4
## 325
         35777_at
                            4
## 326
                            4
         35793_at
```

```
## 327
         35827_at
## 328
         35837_at
                            4
## 329
         35845_at
                            4
## 330 35871_s_at
                            4
## 331
         35877_at
                           13
## 332
         35904_at
                           13
## 333 35939_s_at
                           13
## 334
         35940_at
                           13
## 335
         35949_at
                           13
## 336
                           13
         35972_at
## 337
         35989_at
                           4
## 338
                           4
         35991_at
## 339
         36012_at
                           13
## 340
         36013_at
                           4
## 341
         36017_at
                           13
## 342
         36021_at
                            4
## 343
         36031_at
                           13
## 344
         36046_at
                            4
## 345
         36047_at
                            4
                            4
## 346
         36065_at
                            4
## 347
         36080_at
## 348
         36143_at
                            4
                            4
## 349
         36157_at
## 350
         36188_at
                           13
## 351
         36194_at
                           4
## 352
         36212_at
                           13
## 353
                            4
         36243_at
## 354 36247_f_at
                            4
## 355
         36269_at
                            4
         36274_at
## 356
                           13
## 357
         36358_at
                            4
## 358
         36363_at
                            4
## 359
         36433_at
                            4
## 360 36434_r_at
                            4
## 361
         36510_at
                           13
## 362
                           13
         36521_at
## 363
         36606_at
                            4
## 364
         36622_at
                            4
## 365
         36627_at
                            4
## 366
                           13
         36659_at
                           4
## 367
         36717_at
## 368
         36788_at
                           13
## 369
         367_at
                           13
## 370
         36814_at
                            4
## 371
         36830_at
                           13
## 372
         36913_at
                            4
## 373
                            4
         36914_at
```

```
## 374
        36915_at
## 375
                           4
        36918_at
## 376
         36939_at
                           4
## 377 36968_s_at
                          13
## 378
         36990_at
                           4
## 379
         37006_at
                           4
## 380
         37019_at
                           4
## 381
        37023_at
                          13
## 382
        37056_at
                           4
## 383
                           4
        37058_at
## 384
         37062_at
                           4
## 385
        37067_at
                          13
## 386
        37079_at
                          13
## 387
         37099_at
                          13
## 388
        37109_at
                          13
## 389
        37154_at
                          13
## 390
        37170_at
                           4
## 391
        37172_at
                          13
## 392
        37173_at
                           4
                           4
## 393
         37187_at
                           4
## 394
         37206_at
## 395
        37219_at
                           4
                           4
## 396
        37223_at
## 397
        37243_at
                           4
## 398
       37244_at
                          13
## 399
        37280_at
                           4
## 400
                           4
        37282_at
## 401 37291_r_at
                           4
## 402
         37303_at
                          13
## 403 37322_s_at
                           4
## 404 37323_r_at
                           4
## 405 37356_r_at
                           4
## 406
        37366_at
                           4
## 407
        37404_at
                           4
## 408
       37416_at
                           4
                           4
## 409
       37472_at
                          13
## 410
        37518_at
## 411
       37520_at
                           4
## 412 37521_s_at
                           4
## 413 37522_r_at
                           4
## 414
       37571_at
                          13
## 415
       37578_at
                           4
## 416
       37593_at
                          13
        37619_at
                           4
## 417
## 418
         37658_at
                          13
## 419 37707_i_at
                           4
## 420 37708_r_at
                           4
```

```
## 421
         37723_at
## 422
         37747_at
                            4
## 423
         37748_at
                            4
## 424
                            4
         37752_at
## 425
         37757_at
                           13
## 426
         37767_at
                            4
## 427
         37840_at
                            4
## 428
         37852_at
                            4
## 429
         37926_at
                           13
## 430
         37930_at
                           13
## 431
         37964_at
                            4
## 432
                            4
         38008_at
## 433
         38016_at
                            4
## 434
         38024_at
                            4
## 435 38025_r_at
                            4
## 436
         38035_at
                           13
## 437
         38065_at
                           4
## 438
        38102_at
                           13
## 439
        38120_at
                           4
                            4
## 440
         38168_at
## 441
         38254_at
                           4
## 442 38304_r_at
                           13
## 443
        38353_at
                           13
## 444
         38375_at
                           13
## 445
                            4
         38438_at
## 446
         38485_at
                            4
## 447 38488_s_at
                            4
## 448
         38489_at
                            4
## 449
                            4
        38587_at
## 450
        38606_at
                            4
## 451
        38615_at
                           13
## 452
        38643_at
                           4
## 453
        38649_at
                           13
## 454
         38714_at
                            4
                            4
## 455
         38715_at
## 456
                            4
         38736_at
                            4
## 457 38751_i_at
## 458 38752_r_at
                            4
## 459
         38767_at
                            4
## 460
                            4
         38768_at
                            4
## 461
         38778_at
## 462
        38821_at
                            4
                            4
## 463
        38825_at
## 464
        38838_at
                            4
## 465
         38854_at
                            4
## 466
         38891_at
                            4
## 467
         38957_at
                           13
```

	468	38972_at	13
	469	38988_at	4
##	470	39028_at	13
##	471	39032_at	13
##	472	39037_at	4
##	473	39056_at	4
##	474	39083_at	4
##	475	39131_at	13
##	476	39132_at	4
##	477	39208_i_at	4
##	478	39209_r_at	4
	479	39224_at	4
##	480	39256_at	13
	481	39257_at	13
	482	39269_at	13
		39295_s_at	4
	484	39297_at	13
	485	39333_at	13
	486		4
		39337_at	
	487		4
	488	39369_at	4
	489	39380_at	4
	490	39382_at	4
	491	39405_at	13
		39469_s_at	13
	493	39475_at	4
	494	39481_at	4
	495	39488_at	13
##	496	39489_g_at	13
##	497	39535_at	4
##	498	39536_at	4
##	499	39554_at	4
##	500		4
	501	39576_at	4
	502		13
	503	39600_at	4
	504	39634_at	4
		39662_s_at	4
	506	39665_at	4
	507	39680_at	4
	508	39690_at	4
	509	39698_at	4
	510	39734_at	4
	511	39746_at	4
	512	39748_at	13
		39758_f_at	13
##	514	39777_at	13

10.00		0.077.5	_
	515	39786_at	4
	516	39847_at	4
	517	39850_at	4
##	518	39851_at	4
##	519	39852_at	13
##	520	39878_at	13
##	521	39897_at	4
##	522	39924_at	13
##	523	39929_at	4
	524	39955_at	13
	525	39960_at	4
	526	39979_at	13
	527	40018_at	13
		40058_s_at	4
		40059_r_at	4
		40060_r_at	4
	531	40067_at	13
	532	40072_at	13
	533	40082_at	4
##	534	400_at	13
##	535	40114_at	4
##	536	40121_at	4
##	537		4
##	538	40180_at	13
		40181_f_at	13
	540	40199_at	4
		40217_s_at	4
	542	40218_at	4
	543	40216_at 40225_at	4
	544	40226_at	4
	545	40272_at	4
	546	40310_at	4
	547		13
	548	40323_at	4
	549		4
##	550	40354_at	4
##	551	40392_at	13
##	552	40404_s_at	13
##	553	40449_at	4
	554		4
	555		4
	556	40473_at	13
	557		4
	558		4
	559		13
		40570_at	
		40576_f_at	4
##	561	40633_at	13

```
## 562
      40681_at
                         13
## 563
                          4
       40697_at
## 564
      40710_at
                          4
## 565 40711_at
                          4
## 566
       40727_at
                          4
## 567
        40746_at
                          4
## 568 40770_f_at
                          4
## 569
       40772_at
                          4
## 570
       40773_at
                          4
## 571
      40818_at
                         4
## 572 40828_at
                         13
## 573
      40839_at
                         13
                         4
## 574
      40853_at
## 575 40880_r_at
                         4
## 576
        40893_at
                         13
## 577
          408_at
                         4
## 578 40908_r_at
                         13
## 579
       40943_at
                         4
## 580 40970_at
                         13
                          4
## 581 40989_at
                          4
## 582
       40990_at
## 583 40991_at
                          4
                          4
## 584 40992_s_at
## 585 40993_r_at
                          4
## 586 41014_s_at
                          4
## 587 41024_f_at
                          4
                          4
## 588 41025_r_at
## 589 41026_f_at
                         4
## 590 41069_at
                         13
## 591
      41071_at
                         4
## 592 41104_at
                         4
## 593
      41118_at
                         13
## 594 41119_f_at
                         13
## 595 41145_at
                         4
                         4
## 596 41148_at
## 597 41182_at
                         13
       41191_at
## 598
                         4
## 599
      41276_at
                         13
## 600
       41277_at
                         13
## 601 41300_s_at
                         13
## 602
       41301_at
                         13
## 603
       41308_at
                         4
## 604 41309_g_at
                         4
        41317_at
## 605
                         13
## 606 41318_g_at
                         13
## 607
        41319_at
                         13
                          4
## 608 41376_i_at
```

```
## 609 41377_f_at
## 610
         41391_at
                            4
## 611
         41392_at
                            4
## 612
       41402_at
                            4
## 613
       41434_at
                            4
## 614
        41436_at
                           13
## 615
         41456_at
                            4
## 616
        41459_at
                           13
## 617
         41470_at
                            4
## 618 41491_s_at
                           13
## 619 41492_r_at
                           13
## 620
         41493_at
                           13
## 621
                            4
         41534_at
## 622
         41555_at
                            4
## 623 41556_s_at
                            4
## 624
         41585_at
                            4
## 625 41667_s_at
                           13
## 626 41668_r_at
                           13
## 627
         41697_at
                            4
                            4
## 628
         41801_at
## 629
                            4
         41806_at
## 630
         41860_at
                           13
## 631
           431_at
                            4
## 632
           504_at
                            4
## 633
                            4
         507_s_at
## 634
           579_at
                            4
## 635
                            4
           618_at
## 636
           630_at
                            4
## 637
                            4
         631_g_at
## 638
           655_at
                            4
## 639
         690_s_at
                            4
## 640
         692_s_at
                            4
## 641
         764_s_at
                            4
## 642
           820_at
                            4
                            4
## 643
           886_at
## 644
           931_at
                           13
## 645
                            4
         936_s_at
## 646
         948_s_at
                            4
## 647
           963_at
                           13
## 648
           975_at
                            4
## 649
           990_at
                           13
## 650
         991_g_at
                           13
```

To get this in the classic named-list format:

```
z <- as.list(revmap(x)[chroms])
names(z)</pre>
```

```
## [1] "4" "13"
z[["Y"]]
## NULL
```

Many of the common methods for accessing *Bimap* objects return things in list format. This can be convenient. But you have to be careful about this if you want to use unlist(). For example the following will return multiple probes for each chromosome:

```
chrs = c("12","6")
mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA)
## $`12`
## [1] "1018_at" "1019_g_at" "101_at" "1021_at"
##
## $`6`
## [1] "1026_s_at" "1027_at"
```

But look what happens here if we try to unlist that:

```
unlist(mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA))
## 121 122 123 124 61 62
## "1018_at" "1019_g_at" "101_at" "1021_at" "1026_s_at" "1027_at"
```

Yuck! One trick that will sometimes help is to use Rfunctionunlist2. But be careful here too. Depending on what step comes next, Rfunctionunlist2 may not really help you...

```
unlist2(mget(chrs, revmap(hgu95av2CHR[1:30]), ifnotfound=NA))
## 12 12 12 12 6 6
## "1018_at" "1019_g_at" "101_at" "1021_at" "1026_s_at" "1027_at"
```

Lets ask if the probes in 'pbids' mapped to cytogenetic location "18q11.2"?

To coerce this map to a named vector:

The coercion of the reverse map works too but issues a warning because of the duplicated names for the

reasons stated above:

```
cyto2pb <- as.character(revmap(x))
## Warning: returned vector has duplicated names</pre>
```

2.0.8 Accessing probes that map to multiple targets

In many probe packages, some probes are known to map to multiple genes. The reasons for this can be biological as happens in the arabidopsis packages, but usually it is due to the fact that the genome builds that chip platforms were based on were less stable than desired. Thus what may have originally been a probe designed to measure one thing can end up measuring many things. Usually you don't want to use probes like this, because if they manufacturer doesn't know what they map to then their usefullness is definitely suspect. For this reason, by default all chip packages will normally hide such probes in the standard mappings. But sometimes you may want access to the answers that the manufacturer says such a probe will map to. In such cases, you will want to use the toggleProbes method. To use this method, just call it on a standard mapping and copy the result into a new mapping (you cannot alter the original mapping). Then treat the new mapping as you would any other mapping.

```
## How many probes?
dim(hgu95av2ENTREZID)

## [1] 11533     2

## Make a mapping with multiple probes exposed
multi <- toggleProbes(hgu95av2ENTREZID, "all")
## How many probes?
dim(multi)

## [1] 13393     2</pre>
```

If you then decide that you want to make a mapping that has only multiple mappings or you wish to revert one of your maps back to the default state of only showing the single mappings then you can use toggleProbes to switch back and forth.

```
## Make a mapping with ONLY multiple probes exposed
multiOnly <- toggleProbes(multi, "multiple")
## How many probes?
dim(multiOnly)

## [1] 1860    2

## Then make a mapping with ONLY single mapping probes
singleOnly <- toggleProbes(multiOnly, "single")
## How many probes?
dim(singleOnly)

## [1] 11533    2</pre>
```

Finally, there are also a pair of test methods hasMultiProbes and hasSingleProbes that can be used to see what methods a mapping presently has exposed.

```
## Test the multiOnly mapping
hasMultiProbes(multiOnly)

## [1] TRUE
hasSingleProbes(multiOnly)

## [1] FALSE
## Test the singleOnly mapping
hasMultiProbes(singleOnly)

## [1] FALSE
hasSingleProbes(singleOnly)
```

2.0.9 Using SQL to access things directly

While the mapping objects provide a lot of convenience, sometimes there are definite benefits to writing a simple SQL query. But in order to do this, it is necessary to know a few things. The 1st thing you will need to know is some SQL. Fortunately, it is quite easy to learn enough basic SQL to get stuff out of a database. Here are 4 basic SQL things that you may find handy:

First, you need to know about SELECT statements. A simple example would look something like this:

SELECT * FROM genes;

Which would select everything from the genes table.

SELECT gene_id FROM genes;

Will select only the gene_id field from the genes table.

Second you need to know about WHERE clauses:

SELECT gene_id,_id FROM genes WHERE gene_id=1;

Will only get records from the genes table where the gene_id is = 1.

Thirdly, you will want to know about an inner join:

SELECT * FROM genes, chromosomes WHERE genes._id=chromosomes._id;

This is only slightly more complicated to understand. Here we want to get all the records that are in both the 'genes' and 'chromosomes' tables, but we only want ones where the '_id' field is identical. This is known as an inner join because we only want the elements that are in both of these tables with respect to '_id'. There are other kinds of joins that are worth learning about, but most of the time, this is all you will need to do.

Finally, it is worthwhile to learn about the AS keyword which is useful for making long queries easier to read. For the previous example, we could have written it this way to save space:

SELECT * FROM genes AS g,chromosomes AS c WHERE g._id=c._id;

In a simple example like this you might not see a lot of savings from using AS, so lets consider what happens when we want to also specify which fields we want:

SELECT g.gene_id,c.chromosome FROM genes AS g,chromosomes AS c WHERE g._id=c._id;

Now you are most of the way there to being able to query the databases directly. The only other thing you need to know is a little bit about how to access these databases from R. With each package, you will also get a method that will print the schema for its database, you can view this to see what sorts of tables are present etc.

```
org.Hs.eg_dbschema()
```

To access the data in a database, you will need to connect to it. Fortunately, each package will automatically give you a connection object to that database when it loads.

```
org.Hs.eg_dbconn()
```

You can use this connection object like this:

```
query <- "SELECT gene_id FROM genes LIMIT 10;"
result = dbGetQuery(org.Hs.eg_dbconn(), query)
result</pre>
```

Exercise 5

Retrieve the entrez gene ID and chromosome by using a database query. Show how you could do the same thing by using toTable

2.0.10 Combining data from multiple annotation packages at the SQL level

For a more complex example, consider the task of obtaining all gene symbols which are probed on a chip that have at least one GO BP ID annotation with evidence code IMP, IGI, IPI, or IDA. Here is one way to extract this using the environment-based packages:

```
## Obtain SYMBOLS with at least one GO BP
## annotation with evidence IMP, IGI, IPI, or IDA.
system.time({
bpids <- eapply(hgu95av2GO, function(x) {</pre>
    if (length(x) == 1 && is.na(x))
      NA
    else {
        sapply(x, function(z) {
             if (z$Ontology == "BP")
               z$GOID
             else
               NA
             })
})
bpids <- unlist(bpids)</pre>
bpids <- unique(bpids[!is.na(bpids)])</pre>
g2p <- mget(bpids, hgu95av2G02PR0BE)</pre>
wantedp <- lapply(g2p, function(x) {</pre>
    x[names(x) %in% c("IMP", "IGI", "IPI", "IDA")]
})
```

```
wantedp <- wantedp[sapply(wantedp, length) > 0]
wantedp <- unique(unlist(wantedp))
ans <- unlist(mget(wantedp, hgu95av2SYMBOL))
})
length(ans)
ans[1:10]</pre>
```

All of the above code could have been reduced to a single SQL query with the SQLite-based packages. But to put together this query, you would need to look 1st at the schema to know what tables are present:

```
hgu95av2_dbschema()
```

This function will give you an output of all the create table statements that were used to generate the hgu95av2 database. In this case, this is a chip package, so you will also need to see the schema for the organism package that it depends on. To learn what package it depends on, look at the ORGPKG value:

```
hgu95av20RGPKG
```

Then you can see that schema by looking at its schema method:

dbGetQuery(hgu95av2_dbconn(), "DETACH orgDB"

```
org.Hs.eg_dbschema()
```

So now we can see that we want to connect the data in the go_bp, and symbol tables from the org.Hs.eg.sqlite database along with the probes data in the hgu95av2.sqlite database. How can we do that?

It turns out that one of the great conveniences of SQLite is that it allows other databases to be 'ATTACHed'. Thus, we can keep our data in many differnt databases, and then 'ATTACH' them to each other in a modular fashion. The databases for a given build have been built together and frozen into a single version specifically to allow this sort of behavoir. To use this feature, the SQLite ATTACH command requires the filename for the database file on your filesystem. Fortunately, R provides a nice system independent way of getting that information. Note that the name of the database is always the same as the name of the package, with the suffix '.sqlite':

```
orgDBLoc = system.file("extdata", "org.Hs.eg.sqlite", package="org.Hs.eg.db")
attachSQL = paste("ATTACH '", orgDBLoc, "' AS orgDB;", sep = "")
dbGetQuery(hgu95av2_dbconn(), attachSQL)
## NULL
```

Finally, you can assemble a cross-db sql query and use the helper function as follows. Note that when we want to refer to tables in the attached database, we have to use the 'orgDB' prefix that we specified in the 'ATTACH' query above.:

```
system.time({
SQL <- "SELECT DISTINCT probe_id,symbol FROM probes, orgDB.gene_info AS gi, orgDB.genes AS g, org
zz <- dbGetQuery(hgu95av2_dbconn(), SQL)
})
## user system elapsed
## 0.493 0.036 0.533
#its a good idea to always DETACH your database when you are finished...</pre>
```

NULL

Exercise 6

Retrieve the entrez gene ID, chromosome location information and cytoband information by using a single database query.

Exercise 7

Expand on the example in the text above to combine data from the hgu95av2.db and org.Hs.eg.db with the GO.db package so as to include the GO ID, and term definition in the output.

The version number of R and packages loaded for generating the vignette were:

```
## R version 3.1.1 (2014-07-10)
## Platform: x86_64-apple-darwin10.8.0 (64-bit)
##
## locale:
## [1] en_US.UTF-8/en_US.UTF-8/en_US.UTF-8/C/en_US.UTF-8/en_US.UTF-8
##
## attached base packages:
## [1] parallel stats
                         graphics grDevices utils
                                                     datasets
## [7] methods
               base
##
## other attached packages:
## [1] GO.db_2.14.0
                            hgu95av2.db_2.14.0 AnnotationForge_1.6.1
## [4] org.Hs.eg.db_2.14.0 RSQLite_0.11.4
                                                 DBI_0.3.1
## [7] AnnotationDbi_1.26.1 GenomeInfoDb_1.0.2 Biobase_2.24.0
## [10] BiocGenerics_0.10.0
                            knitr_1.6
##
## loaded via a namespace (and not attached):
## [1] BiocStyle_1.2.0 evaluate_0.5.5 formatR_1.0
                                                    highr_0.3
## [5] IRanges_1.22.10 stats4_3.1.1 stringr_0.6.2 tools_3.1.1
```